

Product datasheet for RC212948L3

SLC6A19 (NM_001003841) Human Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SLC6A19 (NM_001003841) Human Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	SLC6A19
Synonyms:	B0AT1; HND
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC212948).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.

ACCN:	NM_001003841
ORF Size:	1902 bp



[View online »](#)

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001003841.1](#), [NP_001003841.1](#)

RefSeq Size: 1905 bp

RefSeq ORF: 1905 bp

Locus ID: 340024

UniProt ID: [Q695T7](#)


Cytogenetics: 5p15.33

Protein Families: Druggable Genome, Transmembrane

MW: 70.9 kDa

Gene Summary:

This gene encodes a system B(0) transmembrane protein that actively transports most neutral amino acids across the apical membrane of epithelial cells. Mutations in this gene may result in Hartnup disorder, an inherited disease with symptoms such as pellagra, cerebellar ataxia, and psychosis. The expression and function of B0AT1 (SLC6A19) in intestinal cells depends on the presence of the accessory protein angiotensin-converting enzyme 2 (ACE2) which, among other functions, acts as a chaperone for membrane trafficking of B0AT1. The ACE2 is also the cellular receptor for severe acute respiratory syndrome-coronavirus (SARS-CoV) and for SARS-CoV-2 that is causing the coronavirus 2019 (COVID-19) pandemic [provided by RefSeq, Jul 2020]

Product images:

Circular map for RC212948L3