

Product datasheet for **RC212953**

AMPD2 (NM_203404) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	AMPD2 (NM_203404) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	AMPD2
Synonyms:	PCH9; SPG63
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC212953 representing NM_203404
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCTGACCTTCCTTCCCTCCCCCAGGAGCTGTTCAACCCTCACTGGCTGAGAGCGAGCTCCGTAGTG
 CCCGTATGAGTTCCCCGAGGAGAGCCCATTGAACAGCTGGAGGAGCGGCGGCAGCGGCTGGAGCGGCA
 GATCAGCCAGGATGTCAAGCTGGAGCCAGACATCCTGCTTCGGGCCAAGCAAGATTTCTGAAGACGGAC
 AGTGACTCGGACCTACAGCTCTACAAGGAACAGGGTGAGGGGCAGGGTGACCGGAGCCTGCGGGAGCGTG
 ATGTGCTGGAACGGGAGTTTCAGCGGGTACCATTCTGTTGGGAGGAGAAGTGTGGGTGCCGTTACAGA
 CCTGCTGGATGCAGCCAAGAGTGTGTGCGGGCGCTTTCATCCGGGAGAAGTACATGGCCCTGTCCCTG
 CAGAGCTTCTGCCCCACCACCCCGCTACCTGCAGCAGCTGGCTGAAAAGCCTCTGGAGACCCGGACCT
 ATGAACAGGGCCCCGACACCCCTGTGTCTGCTGATGCCCCGGTGACCCCTGCGCTGGAGCAGCACCC
 GTATGAGCACTGTGAGCCAAGCACCATTGCTGGGACCTGGGCTTGGGTCTGCGCATGGTGGGGGTGTG
 GTGCACGTCTACACCCGAGGGAACCCGACGAGCATTGCTCAGAGGTGGAGCTGCCATACCTGACCTGC
 AGGAATTTGTGGCTGACGTCAATGTGCTGATGGCCCTGATTATCAATGGCCCCATAAAGTCACTTCTGCTA
 CCGCCGGCTGCAGTACCTGAGCTCCAAGTTCAGATGCATGTGCTACTCAATGAGATGAAGGAGCTGGCC
 GCCCAGAAGAAAGTGCCACACCGAGATTTCTACAACATCCGCAAGGTGGACCCACATCCATGCCTCGT
 CCTGCATGAACAGAAGCATCTGCTGCGCTTCAAGCGGGCAATGAAGCGGCACCTGGAGGAGATCGT
 GCACGTGGAGCAGGGCCGTGAACAGACGCTGCGGGAGGTCTTTGAGAGCATGAATCTCACGGCTACGAC
 CTGAGTGTGGACACGCTGGATGTGCATGCGGACAGGAACACTTTCATCGCTTTGACAAGTTAATGCCA
 AATAACAACCTATTGGGGAGTCCGCTCCGAGAGATCTTCATCAAGACGGAACAACAGGGTATCTGGGAA
 GTACTTTGCTCACATCATCAAGGAGGTGATGTCAGACCTGGAGGAGAGCAAATACCAGAATGCAGAGCTG
 CGGCTCTCCATTTACGGGCGCTCGAGGGATGAGTGGGACAAGCTGGCGCGCTGGGCCGCTATGCACCGCG
 TGCACTCCCCAACGTGCGCTGGCTGGTGCAGGTGCCCGCCTCTTTGATGTGTACCGTACCAAGGGCCA
 GCTGGCCAACTTCCAGGAGATGCTGGAGAACATCTTCTGCCACTGTTGAGGCCACTGTGCACCCTGCC
 AGCCACCCGAACTGCATCTTCTTAGAGCACGTGGATGGTTTTGACAGCGTGGATGATGAGTCCAAGC
 CTGAAAACCATGTCTTCAACCTGGAGAGCCCCCTGCCTGAGGCGTGGGTGGAGGAGACAACCCACCCTA
 TGCTACTACCTGTACTACCTTTGCCAATGGCCATGTTGAACCACCTGCGCAGGCAGAGGGGCTTC
 CACACGTTTGTGCTGAGGCCACACTGTGGGGAGGCTGGGCCATCCACCCTGGTGTGAGCCTTATGC
 TGGCTGAGAACATTTCCACGGGCTCCTTCTGCGCAAGGCCCGCTCCTGCACTGACTGTACTACCTGGC
 CCAGATCGGCATCGCCATGTCTCCGCTCAGCAACAACAGCCTCTTCTCAGCTATCACCGGAATCCGCTA
 CCGGAGTACCTGTCCCGCGGCTCATGGTCTCCCTGTCCACTGATGATCCCTTGCAGTTCACCTTACCA
 AGGAGCCGCTGATGGAGGAGTACAGCATCGCCACCCAGGTGTGGAAGCTCAGCTCCTGCGATATGTGTGA
 GCTGGCCCGCAACAGCGTGTCTATGAGCGGCTTCTCGACAAGGTAAGAGCCACTGGCTGGGACCCAAC
 TATACCAAGGAAGGCCCTGAGGGGAATGACATCCGCCGACCAATGTGCCAGACATCCGCGTGGGCTACC
 GCTACGAGACCTGTGCCAGGAGCTGGCGCTCATCACGAGGCAGTCCAGAGTGAATGCTGGAGACCAT
 TCCAGAGGAGCGGGTATCACCATGAGCCAGGGCCTCAA

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC212953 representing NM_203404
 Red=Cloning site Green=Tags(s)


MLTFLPSPQELFTRSLAESELSAPYEFPEESPIEQLEERRQRLERQISQDVKLEPDILLRAKQDFLKT
 SDSDLQLYKEQEGQGDRLRERDVLEREFQRTISGEEKCGVPFTDLLDAAKSVVRLFIREKYMALSL
 QSF CPTTRRYLQQLAEKPLETRTYEQGPDTPVSADAPVHPPALEQHPYEHCEPSTMPGDLGLGRMVRGV
 VHVVYTRREPDEHCSEVELPYPDLQEFVADVNVLMALIINGPIKSF CYRRLQYLSSKFQMHVLLNEMKELA
 AQKKVPHRDFYNIKRVKVDTHIHASSCMNQKHLRFIKRAMKRHLEEIVHVEQGREQTLREVFE SMNLTAYD
 LSVDTLDVHADRNTHFRFDKFNAYNP IGESVLR E I F I K T D N R V S G K Y F A H I I K E V M S D L E E S K Y Q N A E L
 RLSIYGRSRDEWDLARWAVMHRVHSPNVRWL VQVPR L F D V Y R T K G Q L A N F Q E M L E N I F L P L F E A T V H P A
 SHPELHLFLEHVDGFSDVDESKPENHVFNLESPLPEAWVEEDNPPYAYLYYTFANMAMLNHLRRQRGF
 HTFVLRPHCGEAGPIHHLVSAFMAENISHGLLLRKAPVLQYL YLAQIGIAMSPLSNNSLFLSYHRNPL
 PEYLSRGLMVSLSTDDPLQHFHTKEPLMEEYSIATQVWKLSSCDMCELARNSVLSMGSFHKVKSHWLGN
 YTKEGPEGNDIRRTNVPDIRVGYRYETLCQELALITQAVQSEMLETIPEEAGITMSPGPQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_203404

ORF Size: 2280 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_203404.1](#), [NP_981949.1](#)

RefSeq Size: 3423 bp

RefSeq ORF: 2283 bp

Locus ID: 271

UniProt ID: [Q01433](#)

Cytogenetics: 1p13.3


Protein Families: Druggable Genome

Protein Pathways: Metabolic pathways, Purine metabolism

MW: 88 kDa

Gene Summary: The protein encoded by this gene is important in purine metabolism by converting AMP to IMP. The encoded protein, which acts as a homotetramer, is one of three AMP deaminases found in mammals. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Apr 2012]

Product images:



Circular map for RC212953