

Product datasheet for **RC212959L3V**

FLRT3 (NM_013281) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	FLRT3 (NM_013281) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FLRT3
Synonyms:	HH21
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_013281
ORF Size:	1947 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC212959).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_013281.2
RefSeq Size:	3773 bp
RefSeq ORF:	1950 bp
Locus ID:	23767
UniProt ID:	Q9NZU0
Cytogenetics:	20p12.1
Domains:	LRRNT, LRRCT, LRR, LRR_TYP, FN3
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)

MW: 72.8 kDa

Gene Summary: This gene encodes a member of the fibronectin leucine rich transmembrane protein (FLRT) family. FLRTs may function in cell adhesion and/or receptor signalling. Their protein structures resemble small leucine-rich proteoglycans found in the extracellular matrix. This gene is expressed in many tissues. Two alternatively spliced transcript variants encoding the same protein have been described for this gene. [provided by RefSeq, Jul 2010]