

## Product datasheet for RC212961

### TLX1NB (NM\_001085398) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** TLX1NB (NM\_001085398) Human Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** TLX1NB  
**Synonyms:** APT-B7; TD1; TDI  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >RC212961 representing NM\_001085398  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCCCGGATCGCC

ATGGCACCCCTCCAATTGCTCCATCCCTGGGGAAGGCGGGGGGCATGGGCAAGTCAGCACTCTCTCCTTT  
 CCCAAGAGGCCATGGGGCTGGGGAGGGGGCAGAGCCTACTCCCTGGGGTACTTGAGGGCAGAGCACTG  
 GGGAGCCTCCCGAGTGGGGCCCCACTTCCCTTTGCCCGGATGAGTGTTCTCTGACACCAGGTTC  
 CCGCCCCACCCTCTAAACAAACAGAAGCATCACCTCCCCACCCTTCTCAAACCCAGAAGGATCTGGTGC  
 CCAGGTCACCTCAGCTGGAGAAATCCAGAATCCGCCTACGGAGAACACTGAGGAACCTTGGTGGAGGCAG  
 GGGTCAGAGGGGACAG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC212961 representing NM\_001085398  
 Red=Cloning site Green=Tags(s)

MAPLQLLHPWRRGAWASQHSLLSQEAMGPGEAEPTPWGYLRAEHWGASPSGAPTLPFAPDECSLTPGS  
 PPHPLNKQKHHPHPSQTQKDLVPRSPQLEKSRIRLRLRNLGGGRGQRGQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mk8023\\_b01.zip](https://cdn.origene.com/chromatograms/mk8023_b01.zip)

**Restriction Sites:** SgfI-MluI



[View online »](#)

**Cloning Scheme:**


**ACCN:** NM\_001085398

**ORF Size:** 366 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001085398.1](#), [NP\\_001078867.1](#)

**RefSeq Size:** 1882 bp


**RefSeq ORF:** 369 bp

**Locus ID:** 100038246

**Cytogenetics:** 10q24.31

MW: 13.2 kDa

**Product images:**



Circular map for RC212961