

## Product datasheet for **RC212971L4V**

### TAF1 48 (TAF1A) (NM\_005681) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	TAF1 48 (TAF1A) (NM_005681) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TAF1 48
Synonyms:	MGC:17061; RAFI48; SL1; TAFI48
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_005681
ORF Size:	1350 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC212971).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_005681.2</a>
RefSeq Size:	1893 bp
RefSeq ORF:	1353 bp
Locus ID:	9015
UniProt ID:	<a href="#">Q15573</a>
Cytogenetics:	1q41
Protein Families:	Transcription Factors
MW:	52.5 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a subunit of the RNA polymerase I complex, Selectivity Factor I (SLI). The encoded protein is a TATA box-binding protein-associated factor that plays a role in the assembly of the RNA polymerase I preinitiation complex. Alternate splicing results in multiple transcript variants encoding multiple isoforms.[provided by RefSeq, Jan 2011]