

Product datasheet for **RC213024A1V**

Human C20orf30 (TMEM230) (NM_001009924) AAV Particle

Product data:

Product Type: AAV Particles

Product Name: Human C20orf30 (TMEM230) (NM_001009924) AAV Particle

Tag: Myc-DDK

Symbol: C20orf30

Synonyms: C20orf30; dj1116H23.2.1; HSPC274

Mammalian Cell Selection: None

Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC213024 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGATGCCGTCCCGTACCAACCTGGCTACTGGAATCCCCAGTAGTAAAGTGAAATATTCAAGGCTCTCCA
GCACAGACGATGGCTACATTGACCTTCAGTTTAAGAAAACCCCTCCTAAGATCCCTTATAAGGCCATCGC
ACTTGCCACTGTGCTGTTTTGATTGGCGCCTTTCTCATTATTATAGGCTCCCTCCTGCTGTCAGGCTAC
ATCAGCAAAGGGGGGCGAGACCGGGCCGTTCCAGTGCTGATCATTGGCATTCTGGTGTTCTACCCGGAT
TTTACCACCTGCGCATCGCTTACTATGCATCCAAGGCTACCGTGCTTACTCCTATGATGACATTCCAGA
CTTTGATGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC213024 protein sequence
Red=Cloning site Green=Tags(s)

MMPSRNLTATGIPSSKVYSRLSSTDGIDYIDLFQFKTPPKIPYKAIALATVLFILGAFLLIIGSLLLSGY
ISKGGADRAVPVLIIGILVFLPGFYHLRIAYYASKGYRAYSYDDIPDFDD

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Species: Human

Serotype: AAV-2

ACCN: NM_001009924



[View online »](#)

ORF Size:	360 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001009924.1, NP_001009924.1</u>
RefSeq Size:	1792 bp
RefSeq ORF:	363 bp
Locus ID:	29058
UniProt ID:	<u>Q96A57</u>
Cytogenetics:	20p13-p12.3
MW:	13.2 kDa