

C20orf30 (TMEM230) (NM_001009924) Human Tagged Lenti ORF Clone

Product Type:	Expression Plasmids
Product Name:	C20orf30 (TMEM230) (NM_001009924) Human Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	C20orf30
Synonyms:	C20orf30; dj11116H23.2.1; HSPC274
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213024).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

... *SgfI* GCG ATC GC C ATG --- // --- NNN⁺ Mlu I ACG CGT ...

ORF

[illegible]


ACCN: NM_001009924

ORF Size: 360 bp

[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001009924.1 , NP_001009924.1
RefSeq Size:	1792 bp
RefSeq ORF:	363 bp
Locus ID:	29058
UniProt ID:	Q96A57
Cytogenetics:	20p13-p12.3
Protein Families:	Transmembrane
MW:	13.2 kDa
Gene Summary:	<p>This gene encodes a multi-pass transmembrane protein that belongs to the TMEM134/TMEM230 protein family. The encoded protein localizes to secretory and recycling vesicle in the neuron and may be involved in synaptic vesicles trafficking and recycling. Mutations in this gene may be linked to familial Parkinson's disease. [provided by RefSeq, Mar 2017]</p>

Product images:



Circular map for RC213024L4