

## Product datasheet for **RC213029L4V**

### ALDH8A1 (NM\_022568) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ALDH8A1 (NM_022568) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ALDH8A1
Synonyms:	ALDH12; DJ352A20.2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_022568
ORF Size:	1461 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213029).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_022568.2</a>
RefSeq Size:	2567 bp
RefSeq ORF:	1464 bp
Locus ID:	64577
UniProt ID:	<a href="#">Q9H2A2</a>
Cytogenetics:	6q23.3
Domains:	aldedh
Protein Families:	Druggable Genome



[View online »](#)

**MW:** 53.4 kDa

**Gene Summary:** This gene encodes a member of the aldehyde dehydrogenase family of proteins. The encoded protein has been implicated in the synthesis of 9-cis-retinoic acid and in the breakdown of the amino acid tryptophan. This enzyme converts 9-cis-retinal into the retinoid X receptor ligand 9-cis-retinoic acid, and has approximately 40-fold higher activity with 9-cis-retinal than with all-trans-retinal. In addition, this enzyme has been shown to catalyze the conversion of 2-aminomuconic semialdehyde to 2-aminomuconate in the kynurenine pathway of tryptophan catabolism. [provided by RefSeq, Jul 2018]