

Product datasheet for **RC213031L4V**

LAD1 (NM_005558) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	LAD1 (NM_005558) Human Tagged ORF Clone Lentiviral Particle
Symbol:	LAD1
Synonyms:	LadA
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_005558
ORF Size:	1551 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213031).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_005558.3
RefSeq Size:	2869 bp
RefSeq ORF:	1554 bp
Locus ID:	3898
UniProt ID:	O00515
Cytogenetics:	1q32.1
Protein Families:	Secreted Protein
MW:	57 kDa



[View online »](#)

Gene Summary:

The protein encoded by this gene may be an anchoring filament that is a component of basement membranes. It may contribute to the stability of the association of the epithelial layers with the underlying mesenchyme. [provided by RefSeq, Jul 2008]