

Product datasheet for **RC213036**

TAF6 (NM_005641) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TAF6 (NM_005641) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	TAF6
Synonyms:	ALYUS; MGC:8964; TAF(II)70; TAF(II)80; TAF2E; TAFII-70; TAFII-80; TAFII70; TAFII80; TAFII85
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RC213036 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)


TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCCGCATCGCC

ATGGCTGAGGAGAAGAAGCTGAAGCTTAGCAACTGTGCTGCCCTCGGAGTCCATGAAGTGGTGGCTG
AATCCATGGGCATCGCCAGATTCAGGAGGAGACCTGCCAGCTGCTAACGGATGAGGTGAGTACCGCAT
CAAAGAGATCGCACAGGATGCCTTGAAGTTCATGCACATGGGGAAGCGGAGAGCTCACCACAGTGAC
ATTGACTACGCCTTGAAGCTAAAGAATGTCGAGCCACTCTATGGCTTCCACGCCAGGAGTTCATTCCTT
TCCGCTTCGCCTCTGGTGGGGCCGGGAGCTTTACTTCTATGAGGAGAAGGAGGTTGATCTGAGCGACAT
CATCAATACCCCTCTGCCCGGGTGGCCCTGGACGTCTGCCTCAAAGCTCATTGGCTGAGCATCGAGGGC
TGCCAGCCAGCTATCCCCGAGAACCAGCCAGCTCCCAAAGAGCAACAGAAGGCTGAAGCCACAGAAC
CCCTGAAGTCAGCAAGCCAGGCCAGGAGAAGACGGACCCTGAAGGGCAAAGGTCAAGGGGCCACCAC
AGCCGACGGCAAAGGGAAAGAGAAGAAGCGCCGCTTGTGGAGGGGGCCCCCTTGGACTGAAGCCC
CGGAGCATCCACGAGTTGTCTGTGGAGCAGCAGCTCTACTACAAGGAGATCACCGAGGCCTGCGTGGGCT
CCTGCGAGGCCAAGAGGGCGGAAGCCCTGCAAAGCATTGCCACGGACCCTGGACTGTATCAGATGCTGCC
ACGGTTCAGTACCTTTATCTCGGAGGGGTCGGTGTGAACGTGGTTTACAGAAACCTGGCCCTACTCATC
TACCTGATGCGTATGGTAAAGCGCTGATGGACAACCCACGCTCTATCTAGAAAAATACGTCCATGAGC
TGATCCAGCTGTGATGACCTGCATCGTGGAGCAGAGTGTGCCTGCGACCAGATGTGGACAATCACTG
GGCACTCCGAGACTTTGCTGCCCGCTGGTGGCCAGATCTGCAAGCATTTTAGCACAACCACTAACAAAC
ATCCAGTCCCGGATCACAAGACCTTACCAAGAGCTGGTGGACGAGAAGACGCCCTGGACGACTCGTT
ATGGCTCCATCGCAGGCTTGGCTGAGCTGGGACACGATGTTATCAAGACTCTGATTCTGCCCGGCTGCA
GCAGGAAGGGGAGCGGATCCGCAGTGTGCTGGACGGCCCTGTGCTGAGCAACATTGACCGGATTGGAGCA
GACCATGTGCAGAGCCTCCTGCTGAAACTGTGCTCCTGTTCTGGCAAAGCTGCGCCACCGCCTGACA
ATCAGGACGCCTATCGGGCAGAATTCGGGTCCTTGGGCCCTCCTCTGCTCCCAGGTGGTCAAGGCTCG
GGCCAGGCTGCTCTGCAGGCTCAGCAGGTCAACAGGACCCTGACCATCACGAGCCCCGGCCACG
CTGACCCTCTCGCAGGCCCCACAGCCTGGCCCTCGCACCCCTGGCTTGTGAAGTTCCTGGCTCCATCG
CACTTCTGTCCAGACTGGTGTCTGCACGAGCGGCTGCCCCACCACAGCCTTCCCTCCTCAACCAA
GTTTATTGTAATGTCATCGTCTCCAGCGCCCCATCCACCAGCAGGTCTGTCCCTCAGCACCTCGGCC
CCCGGCTCAGGTTCCACCACCACTTCGCCGCTCACCACCCTGCCAGCGTGCAGCCATCGTCAAGT
TGGTCTCCACCGCCACCACCGCACCCCCAGCACTGCTCCCTCTGGTCTGGGAGTGTCCAGAAGTACAT
CGTGGTCTCACTTCCCCAACAGGGGAGGGCAAAGGAGGCCCCACCTCCCATCCTTCTCCAGTTCCTCCC
CCGGCATCGTCCCGTCCCCACTCAGCGGAGTGCCTTTGTGGGGGAAGCAGGAGGCTGGGGACAGTC
CCCCTCCAGTCCAGGGACTCCAAAAGCCAATGGCTCCAGCCAACTCCGGTCCCCTCAGCCTGCTCC
G


ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_005641.4
RefSeq Size:	2787 bp
RefSeq ORF:	2034 bp
Locus ID:	6878
UniProt ID:	P49848
Cytogenetics:	7q22.1
Domains:	TAF
Protein Families:	Transcription Factors
Protein Pathways:	Basal transcription factors
MW:	72.7 kDa
Gene Summary:	<p>Initiation of transcription by RNA polymerase II requires the activities of more than 70 polypeptides. The protein that coordinates these activities is transcription factor IID (TFIID), which binds to the core promoter to position the polymerase properly, serves as the scaffold for assembly of the remainder of the transcription complex, and acts as a channel for regulatory signals. TFIID is composed of the TATA-binding protein (TBP) and a group of evolutionarily conserved proteins known as TBP-associated factors or TAFs. TAFs may participate in basal transcription, serve as coactivators, function in promoter recognition or modify general transcription factors (GTFs) to facilitate complex assembly and transcription initiation. This gene encodes one of the smaller subunits of TFIID that binds weakly to TBP but strongly to TAF1, the largest subunit of TFIID. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jun 2010]</p>

Product images:



Circular map for RC213036



Western blot validation of overexpression lysate (Cat# [LY417172]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC213036 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).