

Product datasheet for **RC213046L2V**

ZNF532 (NM_018181) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ZNF532 (NM_018181) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZNF532
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_018181
ORF Size:	3903 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213046).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_018181.4
RefSeq Size:	6493 bp
RefSeq ORF:	3906 bp
Locus ID:	55205
UniProt ID:	Q9HCE3
Cytogenetics:	18q21.32
Domains:	zf-C2H2
MW:	141.5 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]


[View online »](#)