

Product datasheet for RC213051A1V

Human Tcl1 (TCL1A) (NM_001098725) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human Tcl1 (TCL1A) (NM_001098725) AAV Particle
Tag: Myc-DDK
Symbol: Tcl1
Synonyms: TCL1
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC213051 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGGCCGAGTGCCCGACACTCGGGGAGGCAGTCACCGACCACCCGACCGCCTGTGGGCCTGGGAGAAGT
 TCGTGTATTTGGACGAGAAGCAGCAGCCTGGCTGCCCTTAACCATCGAGATAAAGGATAGGTTACAGTT
 ACGGGTGCTCTTGGCTCGGGAAGACGTCGTCTGGGGAGGCCTATGACCCCCACCCAGATAGGCCCAAGC
 CTGCTGCCTATCATGTGGCAGCTCTACCCTGATGGACGATACCGATCCTCAGACTCCAGTTTCTGGCGCT
 TAGTGTACCACATCAAGATTGACGGCGTGGAGGACATGCTTCTCGAGCTGCTGCCAGATGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC213051 protein sequence
 Red=Cloning site Green=Tags(s)
 MAECPTLGEAVTDHPDRLWAWKFFVYLDEKQHAWLPLTIEIKDRLQLRVLLRREDVVLGRPMTPTQIGPS
 LLPIMWQLYPDGRYSSDSSFWRLLVYHIKIDGVEDMLLELLPDD

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_001098725
ORF Size: 342 bp


[View online »](#)

Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_001098725.1
RefSeq Size:	1405 bp
RefSeq ORF:	345 bp
Locus ID:	8115
UniProt ID:	P56279
Cytogenetics:	14q32.13
MW:	13.5 kDa