

Product datasheet for **RC213053L4V**

ZNF160 (NM_033288) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ZNF160 (NM_033288) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZNF160
Synonyms:	F11; HKr18; HZF5; KR18
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_033288
ORF Size:	2454 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213053).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_033288.2
RefSeq Size:	4895 bp
RefSeq ORF:	2457 bp
Locus ID:	90338
UniProt ID:	Q9HCG1
Cytogenetics:	19q13.41-q13.42
Domains:	KRAB, zf-C2H2
Protein Families:	Transcription Factors



[View online »](#)

MW: 93.9 kDa

Gene Summary: The protein encoded by this gene is a Kruppel-related zinc finger protein which is characterized by the presence of an N-terminal repressor domain, the Kruppel-associated box (KRAB). The KRAB domain is a potent repressor of transcription; thus this protein may function in transcription regulation. Multiple transcript variants have been found for this gene. [provided by RefSeq, Apr 2016]