

Product datasheet for **RC213064L4V**

DNMT3A (NM_022552) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	DNMT3A (NM_022552) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DNMT3A
Synonyms:	DNMT3A2; HESJAS; M.HsaIIIA; TBRS
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_022552
ORF Size:	2736 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213064).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_022552.3
RefSeq Size:	4324 bp
RefSeq ORF:	2739 bp
Locus ID:	1788
UniProt ID:	Q9Y6K1
Cytogenetics:	2p23.3
Protein Families:	Druggable Genome
Protein Pathways:	Cysteine and methionine metabolism, Metabolic pathways




[View online »](#)

MW: 101.9 kDa

Gene Summary: CpG methylation is an epigenetic modification that is important for embryonic development, imprinting, and X-chromosome inactivation. Studies in mice have demonstrated that DNA methylation is required for mammalian development. This gene encodes a DNA methyltransferase that is thought to function in de novo methylation, rather than maintenance methylation. The protein localizes to the cytoplasm and nucleus and its expression is developmentally regulated. [provided by RefSeq, Mar 2016]

Product images:



[RC213064L4] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC213064L4V particle to overexpress human DNMT3A-mGFP fusion protein.