

Product datasheet for RC213065L3

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Dopamine D2 Receptor (DRD2) (NM_016574) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: Dopamine D2 Receptor (DRD2) (NM_016574) Human Tagged Lenti ORF Clone

Tag: Myc-DDK

Symbol: Dopamine D2 Receptor

Synonyms: D2DR; D2R Mammalian Cell Puromycin

Selection:

Vector: pLenti-C-Myc-DDK-P2A-Puro (PS100092)

E. coli Selection: Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC213065).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





st The last codon before the Stop codon of the ORF.

ACCN: NM_016574

ORF Size: 1242 bp



Dopamine D2 Receptor (DRD2) (NM_016574) Human Tagged Lenti ORF Clone - RC213065L3

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 016574.2</u>

 RefSeq Size:
 2556 bp

 RefSeq ORF:
 1245 bp

 Locus ID:
 1813

 UniProt ID:
 P14416

Cytogenetics: 11q23.2

Protein Families: Druggable Genome, GPCR, Transmembrane

Protein Pathways: Gap junction, Neuroactive ligand-receptor interaction

MW: 47.2 kDa


Gene Summary: This gene encodes the D2 subtype of the dopamine receptor. This G-protein coupled receptor

inhibits adenylyl cyclase activity. A missense mutation in this gene causes myoclonus dystonia; other mutations have been associated with schizophrenia. Alternative splicing of this gene results in two transcript variants encoding different isoforms. A third variant has been described, but it has not been determined whether this form is normal or due to


aberrant splicing. [provided by RefSeq, Jul 2008]



Product images:



Circular map for RC213065L3



Double digestion of RC213065L3 using Sgfl and Mlul $\,$