

## Product datasheet for **RC213080L2V**

### CDKN3 (NM\_005192) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	CDKN3 (NM_005192) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CDKN3
Synonyms:	CDI1; CIP2; KAP; KAP1
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_005192
ORF Size:	636 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213080).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_005192.2</a>
RefSeq Size:	836 bp
RefSeq ORF:	639 bp
Locus ID:	1033
UniProt ID:	<a href="#">Q16667</a>
Cytogenetics:	14q22.2
Protein Families:	Druggable Genome, Phosphatase
MW:	23.6 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene belongs to the dual specificity protein phosphatase family. It was identified as a cyclin-dependent kinase inhibitor, and has been shown to interact with, and dephosphorylate CDK2 kinase, thus prevent the activation of CDK2 kinase. This gene was reported to be deleted, mutated, or overexpressed in several kinds of cancers. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2008]