

## Product datasheet for **RC213083L4V**

### **p18 INK4c (CDKN2C) (NM\_001262) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	p18 INK4c (CDKN2C) (NM_001262) Human Tagged ORF Clone Lentiviral Particle
Symbol:	p18 INK4c
Synonyms:	INK4C; p18; p18-INK4C
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001262
ORF Size:	504 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213083).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001262.2</a> , <a href="#">NP_001253.1</a>
RefSeq Size:	2104 bp
RefSeq ORF:	507 bp
Locus ID:	1031
UniProt ID:	<a href="#">P42773</a>
Cytogenetics:	1p32.3
Protein Families:	Druggable Genome
Protein Pathways:	Cell cycle



[View online »](#)

**MW:** 17.9 kDa

**Gene Summary:** The protein encoded by this gene is a member of the INK4 family of cyclin-dependent kinase inhibitors. This protein has been shown to interact with CDK4 or CDK6, and prevent the activation of the CDK kinases, thus function as a cell growth regulator that controls cell cycle G1 progression. Ectopic expression of this gene was shown to suppress the growth of human cells in a manner that appears to correlate with the presence of a wild-type RB1 function. Studies in the knockout mice suggested the roles of this gene in regulating spermatogenesis, as well as in suppressing tumorigenesis. Two alternatively spliced transcript variants of this gene, which encode an identical protein, have been reported. [provided by RefSeq, Jul 2008]