

## Product datasheet for **RC213084L1V**

### PGAP2 (NM\_014489) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	PGAP2 (NM_014489) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PGAP2
Synonyms:	CWH43-N; FRAG1; HPMRS3; MRT17; MRT21
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_014489
ORF Size:	945 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213084).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_014489.1</a>
RefSeq Size:	2040 bp
RefSeq ORF:	948 bp
Locus ID:	27315
UniProt ID:	<a href="#">Q9UHJ9</a>
Cytogenetics:	11p15.4
Protein Families:	Druggable Genome, Transmembrane
MW:	36.2 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene plays a role in the maturation of glycosylphosphatidylinositol (GPI) anchors on GPI-anchored proteins. Mutations in this gene are associated with an autosomal recessive syndrome characterized by hyperphosphatasia and intellectual disability. [provided by RefSeq, Jul 2017]