

Product datasheet for **RC213160L3V**

RAB6C (NM_032144) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	RAB6C (NM_032144) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RAB6C
Synonyms:	WTH3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_032144
ORF Size:	762 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213160).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_032144.1 , NP_115520.1
RefSeq Size:	3077 bp
RefSeq ORF:	765 bp
Locus ID:	84084
UniProt ID:	Q9H0N0
Cytogenetics:	2q21.1
Domains:	ras, RAN, RAS, RHO, RAB
Protein Families:	Druggable Genome


[View online »](#)

MW: 28.4 kDa

Gene Summary: May be involved in the regulation of centrosome duplication and cell cycle progression.
[UniProtKB/Swiss-Prot Function]