

## Product datasheet for **RC213161L3V**

### ARSA (NM\_001085427) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ARSA (NM_001085427) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ARSA
Synonyms:	ASA; MLD
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001085427
ORF Size:	1521 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213161).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001085427.1</a>
RefSeq Size:	4127 bp
RefSeq ORF:	1530 bp
Locus ID:	410
UniProt ID:	<a href="#">P15289</a>
Cytogenetics:	22q13.33
Protein Families:	Druggable Genome
Protein Pathways:	Lysosome, Sphingolipid metabolism


[View online »](#)

**MW:** 53.6 kDa

**Gene Summary:** The protein encoded by this gene hydrolyzes cerebroside sulfate to cerebroside and sulfate. Defects in this gene lead to metachromatic leucodystrophy (MLD), a progressive demyelination disease which results in a variety of neurological symptoms and ultimately death. Alternatively spliced transcript variants have been described for this gene. [provided by RefSeq, Dec 2010]