

Product datasheet for **RC213203A1V**

Human Bcl2 Binding component 3 (BBC3) (NM_014417) AAV Particle

Product data:

Product Type: AAV Particles

Product Name: Human Bcl2 Binding component 3 (BBC3) (NM_014417) AAV Particle

Tag: Myc-DDK

Symbol: Bcl2 Binding component 3

Synonyms: JFY-1; JFY1; PUMA

Mammalian Cell Selection: None

Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC213203 representing NM_014417
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCCCGCGCACGCCAGGAGGGCAGCTCCCCGAGCCCGTAGAGGGCCTGGCCCGCGACGGCCCGCGCC
CCTTCCCGCTCGGCCGCTGGTGCCCTCGGCAGTGTCTGCGGCCCTGCGAGCCCGGCTGGTGCCGC
CCCCGCGCCCCACCTGCTGCCCCTGCCTACCTCTGCGCCCCACCGCCCCACCGCCGTACCGCC
GCCCTGGGGGTTCCCGCTGGCTGGGGTCCCCGAGCCGGCCCCGAGGCCCGCGCCCGGACGGTCTCT
AGCCCTCGCTCTCGCTGGCGGAGCAGCACCTGGAGTCGCCGTGCCAGCGCCCCGGGGCTCTGGCGGG
CGGTCCCAACCGCGGCCCGGGAGTCCGCGGGGAGGAGGAACAGTGGGCCCGGGAGATCGGGGCCAG
CTGCGGCGGATGGCGGACGACCTCAACGCACAGTACGAGCGGCGGAGACAAGAGGAGCAGCAGCGGCACC
GCCCCTACCCCTGGAGGGTCTGTACAATCTCATCATGGGACTCTGCCCTTACCCAGGGGCCACAGAGC
CCCCGAGATGGAGCCCAAT

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCTGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC213203 representing NM_014417
Red=Cloning site Green=Tags(s)

MARARQEGSSPEPVEGLARDGPRPFPLGRLVPSAVSCGLCEPGLAAAPAPTLLPAAYLCAPTAPPAVTA
ALGGSRWPGGPRSRPRPDGPQPSLSLAEQHLESPVPSAPGALAGGPTQAAPGVRGEEEQWAREIGAQ
LRRMADDLNAQYERRRQEEQQRHRPSPWRVLYNLIMGLLPLPRGHRAPEMEPN*

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human



[View online »](#)

Serotype:	AAV-2
ACCN:	NM_014417
ORF Size:	579 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_014417.2</u>
RefSeq Size:	1840 bp
RefSeq ORF:	582 bp
Locus ID:	27113
UniProt ID:	<u>Q9BXH1</u>
Cytogenetics:	19q13.32
MW:	20.4 kDa