

Product datasheet for **RC213226A1V**

Human KRTAP19-1 (NM_181607) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human KRTAP19-1 (NM_181607) AAV Particle
Tag:	Myc-DDK
Symbol:	KRTAP19-1
Synonyms:	KAP19.1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC213226 representing NM_181607 Red=Cloning site Blue=ORF Green=Tags(s) TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC GCCCGCATCGCC ATGAGTCACTACGGCAGCTACTACGGAGGCCTGGGCTACAGCTGTGGAGGCTTCGGTGGCCTGGGCTATG GCTATGGCTGTGGATGTGGCAGCTTCTGCAGACGGGTTCTGGCTGTGGCTATGGAGGCTACGGATATGG CTCTGGCTTTGGAAGCTACGGATATGGCTCTGGCTTTGGAGGCTACGGATATGGCTCTGGCTTTGGAGGC TATGGATATGGCTGCTGCCGCCCATCGTACAATGGAGGATACGGATTCTCTGGCTTTTAT ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT ACAAGGATGACGACGATAAGGTTTAA
Protein Sequence:	>RC213226 representing NM_181607 Red=Cloning site Green=Tags(s) MSHYGSYYGGLGYSCGFGGLGYGYGCGGSFCRRSGCGYGGYGYGSGFGSYGYGSGFGGYGYGSGFGG YGYGCCRPSYNGGYGFSGFY TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_181607
ORF Size:	270 bp
Buffer:	PBS with 0.001% Pluronic F68



View online »

Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_181607.1</u>
RefSeq Size:	273 bp
RefSeq ORF:	273 bp
Locus ID:	337882
UniProt ID:	<u>Q8IUB9</u>
Cytogenetics:	21q22.11
MW:	8.8 kDa