

Product datasheet for **RC213228**

DATAF1 (DIDO1) (NM_080796) Human Tagged ORF Clone

Product data:

| | |
|---------------------------|------------------------------------------------------------------------|
| Product Type: | Expression Plasmids |
| Product Name: | DATAF1 (DIDO1) (NM_080796) Human Tagged ORF Clone |
| Tag: | Myc-DDK |
| Symbol: | DATAF1 |
| Synonyms: | BYE1; C20orf158; DATAF-1; DATAF1; DIDO2; DIDO3; DIO-1; DIO1; dj885L7.8 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |



[View online »](#)

ORF Nucleotide Sequence:

>RC213228 representing NM_080796
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGACGACAAAGGCGACCCGAGCAATGAGGAGGCACCTAAGGCCATCAAACCCACCAGCAAAGAGTTCA
 GGAAAACATGGGGTTTTCGAAGGACCACTATCGCCAAGCGAGAGGGCGCAGGGGACGCGGAGGCTGACCC
 ACTGGAGCCGCCACCCACAGCAGCAGCTGGGCTGTCCCTGCGGCGCAGTGGGAGGCAGCCAAAGCGC
 ACTGAGCGCTGGAGCAGTTCCTGACCATTGCGCGGCGCCGCGCAGGAGGAGCATGCCTGTCTCCCTGG
 AGGATTCTGGTGAGCCCACGTCTGCCCGCCACAGCGCGAGACAGCCTCCGAGGGCAGCGTGGAAAAG
 CGCTTCTGAGACCAGAAGCGGCCCCAGTCTGCTTCCACAGCTGTGAAGGAACGACCAGCCTCTTCTGAA
 AAGGTGAAAGGAGGGGATGACCACGATGACACCTCCGATAGTGACAGCGATGGCCTGACCTTGAAGAGC
 TTCAGAATCGCCTTCGCAGGAAGCGGGAACAGGAGCCCACTGAGAGGCCCTGAAAGGGATCCAGAGTCG
 CCTGCGGAAGAAGCGCCGGGAGGAGGTCCCGCGAGACTGTGGGCTCCGAGGCCAGTGACTGTGGAG
 GCGCTCTGCCAGTAAGCAGGAGCCCGAGAACGATCAGGGGTTGTGTCCAGGCTGGGAAAGATGACA
 GAGAGAGTAAGTTGGAGGAAAGGCGGCTCAGGACATCAAAGATGAGGAGCCTGGAGACTTGGGCCGACC
 GAAGCCTGAATGTGAGGGTTACGACCCCAACGCCCTGTATTGCATTTGCCCGCAGCCTCACAAACAGG
 TTTATGATTTGCTGTGACCGCTGTGAAGAATGGTTTCAATGGCGATTGTGTGGCATTCTTCTGAGGCTCGAG
 GGAGGCTTTTGGAAAGGAATGGGGAAGACTATACTGCCCAAACGCACCTTCTGCAAGTGCAGGATGA
 GACTCATTAGAAACGGCAGATCAGCAGGAAGCTAAATGGAGACCTGGAGATGCTGATGGCACCGATTGT
 ACAAGTATAGGAACAATAGAGCAGAAGTCTAGCGAAGACCAAGGGATAAAGGGTAGAATGAGAAAGCTG
 CAAATCCAAGTGGCAAGAAGAACTCAAGATCTCCAGCCTGTGATAGAGGCGCCTGGTGCCTCAAAATG
 TATTGGCCCCGGGTGCTGTACGTGGCGAGCCGACTCGGTGTACTGCAGTAATGACTGTATCCTCAAA
 CACGCCGACGACAATGAAGTTTCTAAGCTCAGGTAAGAACAGAAGCCAAAGCCTAAAGAAAAGATGA
 AGATGAAGCCAGAGAAGCCAGTCTTCCGAAATGCGGTGCTCAGGCAGGTATTAATCTCTTCTGTGCA
 CAAGAGACCAGCTCCAGAAAAAAGAGACCACAGTGAAGAAGGCAGTGGTGGTCCCTGCGCGGAGTGAA
 GCACTCGGAAGGAAGCAGCTTGTGAGAGCAGCACGCCGTCGTGGGCGAGCGATCACAATACAATGCAG
 TAAAGCCAGAAAAGACTGCTGCTCCCTCGCCGTCAGTGTGTATAAATGTATGTATCACCTAGGGGTTGG
 CCTCTGGACCCCTCCCTTCTTCTGGATAGCCATCCCTGGGCTGTCCAGGACTGGGAGTTGCAGCT
 TTGTGT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC213228 representing NM_080796
 Red=Cloning site Green=Tags(s)

MDDKGDPSNEEAPKAIKPTSKEFRKTWGFRRRTIAKREGAGDAEADPLEPPPPQQQLGLSLRRSGRQPKR
 TERVEQFLTIARRRRRRSMPVSLSDSGLTLKELQNLRKREQEPTERPLKGIQSRLRKKRREEGPAETVGSEASDVE
 KVKGGDDHDDTSDSDSGLTLKELQNLRKREQEPTERPLKGIQSRLRKKRREEGPAETVGSEASDVE
 GVLPSKQEPENDQGVVSQAGKDDRESKLEKAAQDIKDEEPGDLGRPKPECEGYDPNALYICRQPHNNR
 FMICCDRCEEFHGDVGISEARGRLLRNGEDYICPNCTILVQVQDETHSETADQEQAKWRPGDADGTD
 TSIGTIEQKSSDQGIKRIEKAANPSGKKLKFQPVIEAPGASKCIGPGCCHVAQPDVSYVCSNDCILK
 HAAATMKFLSSGKEQKPKPEKMKMKPEKPSLPKCGAQAGIKISSVHKRPAPEKKETTAVVVPARSE
 ALGKEAACESSTPSWASDHNYNAVKPEKTAAPSPSLLYKCMYHLGVLLDPSRSFWIAIPWACPLGVAA
 LC

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms:

https://cdn.origene.com/chromatograms/mk6811_d10.zip

Restriction Sites:

Sgfl-Mlul

Cloning Scheme:


ACCN: NM_080796

ORF Size: 1686 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_080796.3](#), [NP_542986.1](#)

RefSeq Size: 2800 bp

RefSeq ORF: 1689 bp

Locus ID: 11083

UniProt ID: [Q9BTC0](#)

Cytogenetics: 20q13.33


Domains: PHD

Protein Families: Druggable Genome, Transcription Factors

MW: 61.4 kDa

Gene Summary: Apoptosis, a major form of cell death, is an efficient mechanism for eliminating unwanted cells and is of central importance for development and homeostasis in metazoan animals. In mice, the death inducer-oblierator-1 gene is upregulated by apoptotic signals and encodes a cytoplasmic protein that translocates to the nucleus upon apoptotic signal activation. When overexpressed, the mouse protein induced apoptosis in cell lines growing in vitro. This gene is similar to the mouse gene and therefore is thought to be involved in apoptosis. Alternatively spliced transcripts have been found for this gene, encoding multiple isoforms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RC213228