

## Product datasheet for **RC213232L4V**

### **ABL2 (NM\_007314) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	ABL2 (NM_007314) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ABL2
Synonyms:	ABLL; ARG
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_007314
ORF Size:	3546 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213232).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_007314.1</a>
RefSeq Size:	3849 bp
RefSeq ORF:	3549 bp
Locus ID:	27
UniProt ID:	<a href="#">P42684</a>
Cytogenetics:	1q25.2
Domains:	ptkase, SH2, TyrKc, SH3, S_TKc
Protein Families:	Druggable Genome, Protein Kinase



[View online »](#)

**Protein Pathways:** ErbB signaling pathway, Viral myocarditis

**MW:** 128.2 kDa

**Gene Summary:** This gene encodes a member of the Abelson family of nonreceptor tyrosine protein kinases. The protein is highly similar to the c-abl oncogene 1 protein, including the tyrosine kinase, SH2 and SH3 domains, and it plays a role in cytoskeletal rearrangements through its C-terminal F-actin- and microtubule-binding sequences. This gene is expressed in both normal and tumor cells, and is involved in translocation with the ets variant 6 gene in leukemia. Multiple alternatively spliced transcript variants encoding different protein isoforms have been found for this gene. [provided by RefSeq, Nov 2009]