

## Product datasheet for **RC213233L4V**

### KLHL32 (NM\_052904) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	KLHL32 (NM_052904) Human Tagged ORF Clone Lentiviral Particle
Symbol:	KLHL32
Synonyms:	BKLDH5; dj21F7.1; KIAA1900; UG0030H05
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_052904
ORF Size:	1860 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213233).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_052904.1</a>
RefSeq Size:	3322 bp
RefSeq ORF:	1863 bp
Locus ID:	114792
UniProt ID:	<a href="#">Q96NJ5</a>
Cytogenetics:	6q16.1
Domains:	BTB, Kelch
MW:	70.2 kDa



[View online »](#)