

## Product datasheet for **RC213250L2V**

### IRF2BP2 (NM\_182972) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	IRF2BP2 (NM_182972) Human Tagged ORF Clone Lentiviral Particle
Symbol:	IRF2BP2
Synonyms:	CVID14
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_182972
ORF Size:	1761 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213250).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_182972.1</a> , <a href="#">NP_892017.2</a>
RefSeq Size:	1900 bp
RefSeq ORF:	1764 bp
Locus ID:	359948
UniProt ID:	<a href="#">Q7Z5L9</a>
Cytogenetics:	1q42.3
Protein Families:	Druggable Genome
MW:	60.8 kDa



[View online »](#)

**Gene Summary:**

This gene encodes an interferon regulatory factor-2 (IRF2) binding protein that interacts with the C-terminal transcriptional repression domain of IRF2. Alternative splicing results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, Jul 2008]