

Product datasheet for **RC213252L3V**

AGAP8 (AGAP4) (NM_133446) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	AGAP8 (AGAP4) (NM_133446) Human Tagged ORF Clone Lentiviral Particle
Symbol:	AGAP8
Synonyms:	AGAP-4; AGAP-8; AGAP8; CTGLF1; CTGLF5; MRIP2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_133446
ORF Size:	1989 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213252).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_133446.1 , NP_597703.2
RefSeq Size:	2420 bp
RefSeq ORF:	1992 bp
Locus ID:	119016
UniProt ID:	Q96P64
Cytogenetics:	10q11.22
Domains:	ArfGap, PH, ANK
MW:	72.9 kDa



[View online »](#)

Gene Summary: Putative GTPase-activating protein.[UniProtKB/Swiss-Prot Function]