

## Product datasheet for **RC213264L4V**

### Inactive ubiquitin carboxyl terminal hydrolase 50 (USP50) (NM\_203494) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Inactive ubiquitin carboxyl terminal hydrolase 50 (USP50) (NM_203494) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Inactive ubiquitin carboxyl terminal hydrolase 50
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_203494
ORF Size:	1002 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213264).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_203494.1</a>
RefSeq Size:	1020 bp
RefSeq ORF:	1005 bp
Locus ID:	373509
UniProt ID:	<a href="#">Q70EL3</a>
Cytogenetics:	15q21.2
Protein Families:	Druggable Genome
MW:	38.2 kDa



[View online »](#)

**Gene Summary:** Has no peptidase activity.[UniProtKB/Swiss-Prot Function]