

Product datasheet for RC213287A1V

Human CLEC5A (NM_013252) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human CLEC5A (NM_013252) AAV Particle
Tag: Myc-DDK
Symbol: CLEC5A
Synonyms: CLECSF5; MDL-1; MDL1
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC213287 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGAACTGGCACATGATCATCTCTGGGCTTATTGTGGTAGTGCTTAAAGTTGTTGGAATGACCTTATTTCT
 TACTTTATTTCCACAGATTTTTAACAAAAGTAACGATGGTTTCACCACCACAGGAGCTATGGAACAGT
 CTCACAGATTTTGGGAGCAGTCCCAAGTCCCAACGGCTTCATTACCACAAGGAGCTATGGAACAGTC
 TGCCCCAAGACTGGGAATTTATCAAGCAAGATGTTTTCTTATCCACTTCTGAATCATCTTGAATG
 AAAGCAGGGACTTTTGCAAAGGAAAAGGATCCACATTGGCAATTGTCAACACGCCAGAGAACTGAAGTT
 TCTTCAGGACATAACTGATGCTGAGAAGTATTTTATTGGCTTAATTTACCATCGTGAAGAGAAAAGGTGG
 CGTTGGATCAACAACCTCTGTGTTCAATGGCAATGTTACCAATCAGAATCAGAATTTCAACTGTGCGACAA
 TTGGCCTAACAAAGACATTTGATGCTGCATCATGTGACATCAGCTACCGCAGGATCTGTGAGAAGAATGC
 CAAA

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC213287 protein sequence
 Red=Cloning site Green=Tags(s)

MNWHMIISGLIVVVLKVVGMTLFLLYFPQIFNKSNDGFTTTRSYGTVSQIFGSSSPSPNGFITTRSYGTV
 CPKDWEFYQARCFLLSTESSWNESRDFCKGKSTLAIVNTPEKLKFLQDITDAEKYFIGLIYHREEKRW
 RWINNSVFNGNVTNQNQNFNCATIGLTKTFDAASCDISYRRICEKNAK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human


[View online »](#)

Serotype:	AAV-2
ACCN:	NM_013252
ORF Size:	564 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_013252.2</u>
RefSeq Size:	3524 bp
RefSeq ORF:	567 bp
Locus ID:	23601
UniProt ID:	<u>Q9NY25</u>
Cytogenetics:	7q34
MW:	21.5 kDa