

Product datasheet for RC213295A1V

Human C17orf87 (SCIMP) (NM_207103) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human C17orf87 (SCIMP) (NM_207103) AAV Particle
Tag: Myc-DDK
Symbol: C17orf87
Synonyms: C17orf87; UNQ5783
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC213295 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGATACTTTACAGTTCAGGATTCCTGCAATGAGCTGGTGGAGGAATAATTTCTGGATCATCTTAG
 CTGTGGCCATCATCGTTGTCTCTGTGGGTCTGGGCCTCATCTGTACTGTGTCTGTAAGTGGCAGCTTAG
 ACGAGGCAAGAAATGGGAAATGGCAAGCCCTGAAACACAAGCAAGTAGATGAAGAAAAGATGTATGAG
 AATGTTCTTAATGAGTCGCCAGTTCAATTACCGCCTCTGCCACCGAGGAATTGGCCTTCTCTAGAAGACT
 CTTCCACAGGAAGCCCAAGTCAGCCGCCCTACATACTACTGGTAAATAAAGTTAAAAATAAGAA
 GACTGTTTCCATCCCAAGCTACATTGAGCCTGAAGATGACTATGACGATGTTGAAATCCCTGCAATACT
 GAAAAAGCATCATTT

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC213295 protein sequence
 Red=Cloning site Green=Tags(s)

MDTFTVQDSTAMSWRRNFWIILAVAIIVSVGLGLILYCVCKWQLRRGKKWEIAKPLKHKQVDEEKMYE
 NVLNESPVQLPPLPPRNWPSLEDSSPQEAPSPATYSLVNKVNKKTIVSIPSYIEPEDDYDDVEIPANT
 EKASF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2


[View online »](#)

ACCN:	NM_207103
ORF Size:	435 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_207103.1</u>
RefSeq Size:	2424 bp
RefSeq ORF:	438 bp
Locus ID:	388325
UniProt ID:	<u>Q6UWF3</u>
Cytogenetics:	17p13.2
MW:	16.6 kDa