

## Product datasheet for **RC213295L2V**

### **C17orf87 (SCIMP) (NM\_207103) Human Tagged ORF Clone Lentiviral Particle**

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	C17orf87 (SCIMP) (NM_207103) Human Tagged ORF Clone Lentiviral Particle
Symbol:	C17orf87
Synonyms:	C17orf87; UNQ5783
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_207103
ORF Size:	435 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213295).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_207103.1</a>
RefSeq Size:	2424 bp
RefSeq ORF:	438 bp
Locus ID:	388325
UniProt ID:	<a href="#">Q6UWF3</a>
Cytogenetics:	17p13.2
Protein Families:	Transmembrane
MW:	16.6 kDa


[View online »](#)

**Gene Summary:**

This gene encodes a transmembrane adaptor protein that is expressed in antigen-presenting cells and is localized in the immunologic synapse. The encoded protein is involved in major histocompatibility complex class II signal transduction and immune synapse formation. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Dec 2012]