

## Product datasheet for **RC213303A1V**

### Human ECSCR (NM\_001077693) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Human ECSCR (NM\_001077693) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** ECSCR  
**Synonyms:** ARIA; ECSM2  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >RC213303 representing NM\_001077693  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGCATCGCC**

ATGGGCACCGCAGGAGCCATGCAGCTGTGCTGGGTGATCCTGGGCTTCTCTGTTCGAGGCCACAAC  
CCCAGCCCACAATGACCCAGACCTCTAGCTCTCAGGGAGGCCTTGGCGGTCTAAGCTGACCACAGAGCC  
AGTTTCTTCCAACCCAGGATACATCCCTTCTCAGAGGCTAACAGGCCAAGCCATCTGTCCAGCACTGGT  
ACCCAGGCGCAGGTGTCCCGAGCAGTGAAGAGACGGAGGCACAAGCAGAGACACATTTCAAAGTGTTC  
CCCCAATTCAACCACCATGAGCCTGAGCATGAGGGAAGATGCGACCATCCTGCCAGCCCCACGTCAGA  
GACTGTGCTCACTGTGGCTGCATTTGGTGTTATCAGCTTCATTGTCATCCTGGTGGTGTGGTGATCATC  
CTAGTTGGTGTGGTCAGCCTGAGGTTCAAGTGTGGAAGAGCAAGGAGTCTGAAGATCCCAGAAACCTG  
GGAGTTTCAGGGCTGTCTGAAAGCTGCTCCACAGCCAATGGAGAGAAAGACAGCATCACCTTATCTCCAT  
GAAGAACATCAACATGAATAATGGCAAACAAAGTCTCTCAGCAGAGAAGGTTCTT

**ACGCGT**ACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC213303 representing NM\_001077693  
Red=Cloning site Green=Tags(s)

MGTAGAMQLCWILGFLLFRGHNSQPTMTQTSSSQGGLGGLSLTTEPVSSNPGYIPSSSEANRPSHLSSTG  
TPGAGVPSSGRDGGTSRDTFQTPPNSTTMSLSMREDATILPSPTSETVLTVAAGVISFIVILVVVVII  
LVGVVSLRFKCRKSKESDPQKPGSSGLSESCSTANGEKDSITLISMKNINMNGKQSLSAEKVL

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Species:** Human



[View online »](#)

|               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|---------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Serotype:     | AAV-2                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| ACCN:         | NM_001077693                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| ORF Size:     | 615 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Buffer:       | PBS with 0.001% Pluronic F68                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Stability:    | AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus. |
| RefSeq:       | <u><a href="#">NM_001077693.1</a></u>                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| RefSeq Size:  | 1030 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| RefSeq ORF:   | 618 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Locus ID:     | 641700                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| UniProt ID:   | <u><a href="#">Q19T08</a></u>                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Cytogenetics: | 5q31.2                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| MW:           | 21.1 kDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |