

## Product datasheet for **RC213315L4V**

### PHF20L1 (NM\_032205) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	PHF20L1 (NM_032205) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PHF20L1
Synonyms:	CGI-72; TDRD20B
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_032205
ORF Size:	1719 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213315).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_032205.3</a> , <a href="#">NP_115581.3</a>
RefSeq Size:	2286 bp
RefSeq ORF:	1721 bp
Locus ID:	51105
Cytogenetics:	8q24.22
Protein Families:	Druggable Genome, Transcription Factors
MW:	63.3 kDa



[View online »](#)