

Product datasheet for **RC213337A1V**

Human CPEB1 (NM_001079534) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human CPEB1 (NM_001079534) AAV Particle
Tag:	Myc-DDK
Symbol:	CPEB1
Synonyms:	CPE-BP1; CPEB; CPEB-1; h-CPEB; hCPEB-1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC213337 ORF sequence

Red=Cloning site Blue=ORF Green=Tags(s)

 TTTTGAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

 ATGCTTTTCCCAACCTCTGCGCAAGAATCTTCCCGTGGCCTCCAGATGCAAATGACTTGTGCCTTGCC
TGCAGTCCCTCAGTCTGACAGGCTGGGACCGACCCTGGAGCACCCAGGACTCAGATTCTCAGCCAGAG
CAGCACACACTCGGTACTGAGCATGCTCCATAACCCACTGGGAAATGTCCTAGGAAAACCCCTTGAGC
TTCCTGCCTCTGGATCCCTTGGGTCTGACTTGGTGGACAAGTTTCCAGCACCCCTCAGTTAGAGGATCAC
GCCTGGACACCCGGCCATCCTGGACTCTCGATCTAGCAGCCCTCTGACTCAGACACCACTGGCTTCAG
CTCTGGATCAGATCATCTCTCAGATTTGATTTCAAGCCTTCGATTTCTCCACCTCTGCCCTTCTGTCT
CTGTCAAGGGGTGGTCCAGAGACCTTTAAAGATGGGGTAGGGTCTCGGATGGACCAAGAGCAAGCTG
CTTTGCTGAGTCACTCCCTCCCAACCACTGCTTCAAAGAGATGGCCAGGAGCTTCTGTGTGCCATC
CTGGGACCTCCTGAAGCTCCCAAGACCCCTTCAGCATAGAGAGAGAGGCCAGGCTGCACCGACAAGCT
GCAGCTGTGAATGAAGCCACCTGTACCTGGAGTGGCCAGCTTCTCCCGGAACCTATAAGAACCCCATCT
ACTCTTGAAGGTGTTTCTAGGAGGTGTTCTTGGGATATTACAGAAGCTGGATTAGTTAACACCTTCCG
TGTTTTTGCTCTTTGAGTGTGGAGTGGCTGGTAAGGATGGCAAGCATCCCCGGTGTCTCCAAAGGG
TATGTGTATCTGGTCTTGAAGTGAAGTCTGTCCGATCCTTGCTTCAAGCTTGTCTCATGACCCGC
TGAGCCAGATGGCTGAGTGAATATTATTTCAAGATGTCAGCCGAAGGATGCGCTGCAAGGAGGTGCA
GGTATCCCTGGGTATTAGCCGACAGTAACCTTGTCCGAGCCCATCTCAGAGGCTTGACCCAGCAGG
ACGGTGTGTCGGTGTCTGCATGGAATGCTAAATGCTGAGGCCCTGGCAGCCATCTTGAACGACCTAT
TTGTGGAGTGGTGTATGCCGGGATTGACACAGATAAGCACAAGTATCCATTGGTTCTGGTGTGTGAC
TTTCAATAACCAACGGAGTTACCTGAAAGCAGTCAGCGCTGCTTTTGTGGAGATCAAAACCAACCAAGTTC
ACAAAGAAGGTTCAAGTTGACCCCTACCTGGAAGATTCTCTGTGTATATCTGCAGTTCTCAGCTGGTC
CTTTCTTGTGAGATCAGGTCTGCTTCAAATACTTCTGCCGAGCTGCTGGCACTGGCGGCACAGCAT
GGAGGGCTGCCACCACAGCCCTGATGCGGAACCAAGAACCGAGATTCCAGC

 A**CGCGT**ACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC213337 protein sequence

Red=Cloning site Green=Tags(s)

 MLFPTSAQESSRGLPDANDLCLGLQSLSLTGWDRPWSTQSDSSAQSSSTHSVL SMLHNPLGNVLGKPPLS
FLPLDPLGSDLVDKFPAPSVRGSRLDTRPILDSRSSPSDSTSGFSSGSDHLSLISLRISPPLPFLS
LSGGGPRDPLKMGVGSMDQEQAALAAVTPSPTSASKRWPGASVWPSWDLLEAPKDPFSIEREARLHRQA
AAVNEATCTWSGQLPPRNYKNPIYSCKVFLGGVPWDITEAGLVNTRFVFGSLSVEWPGDKGHPRCPPKG
YVYLVFELEKSVRSLLQACSHDPLSPDGLSEYFYMSSRRMRCKEVQVIPWVLADSNFVRSPSQRLDPSR
TVFVGALHGMLNAEALAILNDLFGGVYAGIDTDKHKYPIGSGRVTFNQRSYLKAVSAFVEIKTTKF
TKKVQIDPYLEDLSLCHICSSQPGPFRCRDQVCFKYFCRSCWHWRHSMEGLRHHSPLMRNQKNRDSS

 TRTRPLE**QKLISEEDLA**NDILDYKDDDDKV

Species:

Human

Serotype:

AAV-2

ACCN:

NM_001079534

ORF Size:

1458 bp

Buffer:

PBS with 0.001% Pluronic F68

Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001079534.1</u> , <u>NP_001073002.1</u>
RefSeq Size:	3095 bp
RefSeq ORF:	1461 bp
Locus ID:	64506
UniProt ID:	<u>Q9BZB8</u>
Cytogenetics:	15q25.2
MW:	53.6 kDa