

## Product datasheet for **RC213341L2V**

### GPD2 (NM\_001083112) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	GPD2 (NM_001083112) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GPD2
Synonyms:	GDH2; GPDM; mGDH; mGPDH
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_001083112
ORF Size:	2181 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213341).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001083112.1</a>
RefSeq Size:	6046 bp
RefSeq ORF:	2184 bp
Locus ID:	2820
UniProt ID:	<a href="#">P43304</a>
Cytogenetics:	2q24.1
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Glycerophospholipid metabolism



[View online »](#)

**MW:** 80.9 kDa

**Gene Summary:** The protein encoded by this gene localizes to the inner mitochondrial membrane and catalyzes the conversion of glycerol-3-phosphate to dihydroxyacetone phosphate, using FAD as a cofactor. Along with GDP1, the encoded protein constitutes the glycerol phosphate shuttle, which reoxidizes NADH formed during glycolysis. Two transcript variants encoding the same protein have been found for this gene.[provided by RefSeq, Jan 2010]