

Product datasheet for **RC213346L4V**

TECTB (NM_058222) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TECTB (NM_058222) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TECTB
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_058222
ORF Size:	987 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213346).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_058222.1
RefSeq Size:	2810 bp
RefSeq ORF:	990 bp
Locus ID:	6975
UniProt ID:	Q96PL2
Cytogenetics:	10q25.2
MW:	37 kDa
Gene Summary:	This gene encodes a non-collagenous glycoprotein component of the tectorial membrane, which covers the auditory hair cells in the cochlea of the inner ear. A similar protein in mouse functions in low-frequency hearing. [provided by RefSeq, Jul 2013]



[View online »](#)