

## Product datasheet for **RC213357A1V**

### Human MRP6 (ABCC6) (NM\_001079528) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Human MRP6 (ABCC6) (NM\_001079528) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** MRP6  
**Synonyms:** ABC34; ARA; EST349056; GACI2; MLP1; MOAT-E; MOATE; MRP6; PXE; PXE1; URG7  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >RC213357 representing NM\_001079528  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGCCGCGCCTGCTGAGCCCTGCGCGGGGCAGGGGTCTGGAACCAGACAGAGCCTGAACCTGCCGCCA  
 CCAGCCTGCTGAGCCTGTGCTTCTGAGAACAGCAGGGGTCTGGGTACCCCCATGTACCTCTGGGTCTT  
 TGGTCCCATCTACCTCCTTTCATCCACCACCATGGCCGGGGCTACCTCCGGATGTCCCACTCTTCAA  
 AAGATGGTAGCTGCCATCCCTGGGAGCCTGGAACCAGGCAATGTTCCGGGGAGGCAGGGGACAGGCT  
 GGAACCTGGTGAAGTCT

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC213357 representing NM\_001079528  
 Red=Cloning site Green=Tags(s)

MAAPAEPGAGQGVWNQTEPEPAATSLLSLCFLRTAGVWVPPMYLWVLGPIYLLFIHHHGRGYLRMSPLFK  
 AKMVAaipgslepgnvrgrqgtgwnlvks

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Species:** Human  
**Serotype:** AAV-2  
**ACCN:** NM\_001079528  
**ORF Size:** 297 bp



[View online »](#)

<b>Buffer:</b>	PBS with 0.001% Pluronic F68
<b>Stability:</b>	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
<b>RefSeq:</b>	<u>NM_001079528.1</u>
<b>RefSeq Size:</b>	727 bp
<b>RefSeq ORF:</b>	300 bp
<b>Locus ID:</b>	368
<b>UniProt ID:</b>	<u>O95255</u>
<b>Cytogenetics:</b>	16p13.11
<b>MW:</b>	10.6 kDa