

Product datasheet for **RC213360L3V**

AKIRIN1 (NM_024595) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	AKIRIN1 (NM_024595) Human Tagged ORF Clone Lentiviral Particle
Symbol:	AKIRIN1
Synonyms:	C1orf108; STRF2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_024595
ORF Size:	576 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213360).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_024595.1
RefSeq Size:	1968 bp
RefSeq ORF:	579 bp
Locus ID:	79647
UniProt ID:	Q9H9L7
Cytogenetics:	1p34.3
MW:	21.7 kDa



[View online »](#)

Gene Summary:

Functions as signal transducer for MSTN during skeletal muscle regeneration and myogenesis. May regulates chemotaxis of both macrophages and myoblasts by reorganising actin cytoskeleton, leading to more efficient lamellipodia formation via a PI3 kinase dependent pathway.[UniProtKB/Swiss-Prot Function]