

Product datasheet for **RC213409L3V**

CTSA (NM_000308) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CTSA (NM_000308) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CTSA
Synonyms:	GLB2; GSL; NGBE; PPCA; PPGB
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_000308
ORF Size:	1491 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213409).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_000308.2 , NP_000299.2
RefSeq Size:	2254 bp
RefSeq ORF:	1443 bp
Locus ID:	5476
UniProt ID:	P10619
Cytogenetics:	20q13.12
Domains:	serine_carbpept
Protein Families:	Druggable Genome, Protease



[View online »](#)

Protein Pathways: Renin-angiotensin system

MW: 56.12 kDa

Gene Summary: This gene encodes a member of the peptidase S10 family of serine carboxypeptidases. Alternative splicing results in multiple transcript variants, at least one of which encodes a preproprotein that is proteolytically processed to generate two chains that comprise the heterodimeric active enzyme. This enzyme possesses deamidase, esterase and carboxypeptidase activities and acts as a scaffold in the lysosomal multienzyme complex. Mutations in this gene are associated with galactosialidosis. [provided by RefSeq, Nov 2015]