

Product datasheet for **RC213422L4V**

Entactin (NID1) (NM_002508) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Entactin (NID1) (NM_002508) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Entactin
Synonyms:	NID
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_002508
ORF Size:	3741 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213422).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_002508.1
RefSeq Size:	4898 bp
RefSeq ORF:	3744 bp
Locus ID:	4811
UniProt ID:	P14543
Cytogenetics:	1q42.3
Domains:	Idl_recept_b, thyroglobulin_1, EGF_CA, NIDO, EGF, EGF, G2F
Protein Families:	Druggable Genome, Secreted Protein



[View online »](#)

MW: 136.4 kDa

Gene Summary: This gene encodes a member of the nidogen family of basement membrane glycoproteins. The protein interacts with several other components of basement membranes, and may play a role in cell interactions with the extracellular matrix. [provided by RefSeq, Jul 2008]