

Product datasheet for RC213426

C22orf9 (KIAA0930) (NM_015264) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C22orf9 (KIAA0930) (NM_015264) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	C22orf9
Synonyms:	C22orf9; LSC3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC213426 representing NM_015264 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGCC

ATGGGCTCTCAGGCTGCTGCTGAGTGGAGGAAGTGGGCTCCTGGGAGGTGTCTCCAGCCTCTCTGGAT
 GCTCCATGGGGTCTTCAAGGATGACCGCATCGTCTTCTGGACTTGGATGTTCTCCACCTACTTCATGGA
 GAAATGGGCTCCCCGGCAGGACGACATGCTTTTCTATGTGCGCCGAAGCTGGCGTACTCCGGCAGCGAA
 AGCGGTGCAGACGGGAGGAAGGCAGCTGAGCCTGAGGTGGAGGTGGAGGTGTACCGCGGGACTCCAAGA
 AGCTGCCAGGCCCTGGGAGACCCTGACATCGACTGGGAGGAGAGCGTCTGCCTGAATCTCATCCTGCAGAA
 GCTGGACTACATGGTGACCTGTGCGGTGTGCACACGTGCTGACGGCGGGGACATTACATCCATAAGAAG
 AAATCTCAGCAAGTGTTGCGCTCCCCAGTAAACACCCCATGGACAGCAAGGGGGAGGAGTCCAAGATCA
 GCTACCCCAACATCTTCTTCATGATTGACAGCTTCGAGGAGGTGTTTACGCGACATGACCGTAGGGGAAGG
 AGAGATGGTCTGTGTGGAGCTGGTGGCTAGTGACAAAACCAACAGTTCAGGGGGTTCATCTTCAGGGC
 TCCATCCGCTACGAGGCGCTCAAGAAGGTGTATGACAACCGGGTGAGCGTGGCCGCCCGCATGGCACAGA
 AGATGTCGTTTGGCTTCTACAAGTACAGCAACATGGAGTTTGTGCGCATGAAGGGCCCCAGGGCAAGGG
 CCACGCCGAGATGGCGGTGAGCGAGTGTCTACAGGTGACACATCCCCCTGTGGGACTGAAGAGGACTCC
 AGCCCAGCTTCGCCCATGCACGAGCGGGTGACCTCCTTCAGCACACCCCCACCCAGAACGGAACAACC
 GGCTGCCTTCTTCTCCCATCCCTCAAGAGGAAGGTGCCCGGAACCGGATCGCTGAGATGAAGAAGTC
 GCACTCGGCCAACGACAGCGAGGAGTCTTCCGGGAGGACGACGGTGGAGCCGATCTGCACAATGCAACC
 AACCTGCGGTCTCGGTCCTGTGCGGCACAGGACGGTCCCTGGTCGGGTCTGGCTGAAGCTGAACAGAG
 CAGATGGAAACTCCTTCTCTATGCACACTTAACCTACGTACGTTGCCGCTGCATCGGATTTTAACAGA
 CATCCTGGAAGTTCGGCAGAAGCCATCCTGATGACC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence: >RC213426 representing NM_015264
 Red=Cloning site Green=Tags(s)

MGSQAAAEWRNWSWEVSSSLSGCSMGCFKDDRIVFWTWMFSTYFMEKWAPRQDDMLFYVRRKLAYSQSE
 SGADGRKAAEPEVEVEVYRRDSKKLPGLGDPDIDWEESVCLNLILQKLDYMTCAVCTRADGGDIHKK
 KSQQVFASPSKHPMDSKGEESKISYPNIFFMIDSFEFVSDMTVGEGEMVCVELVASDKTNTFQGVIFQG
 SIRYEALKKVYDNRVSAARMAQKMSFGFYKYSNMEFVRMKGPQKGHAEMAVSRVSTGDTSPCGTEEDS
 SPASPMHERVTSFSTPPTPERNNRPAFFSPSLKRVKVPNRNIAEMKKSHSANDSEFFREDDGGADLHNAT
 NLRSRSLSGTGRSLVGSWLKLNRADGNFLLYLAHTYVTLPLHRILTDILEVRQKPIILMT

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8023_h06.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_015264

ORF Size: 1227 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: NM_015264.2

RefSeq Size: 6287 bp

RefSeq ORF: 1230 bp


Locus ID: 23313

UniProt ID: Q6ICG6


Cytogenetics: 22q13.31

MW: 46 kDa


Product images:



Circular map for RC213426



Western blot validation of overexpression lysate (Cat# [LY414655]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC213426 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified KIAA0930 protein (Cat# [TP313426]). The protein was produced from HEK293T cells transfected with KIAA0930 cDNA clone (Cat# RC213426) using MegaTran 2.0 (Cat# [TT210002]).