

Product datasheet for **RC213430L3V**

CLUAP1 (NM_024793) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CLUAP1 (NM_024793) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CLUAP1
Synonyms:	CFAP22; FAP22; IFT38
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_024793
ORF Size:	1242 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213430).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_024793.1
RefSeq Size:	3616 bp
RefSeq ORF:	744 bp
Locus ID:	23059
UniProt ID:	Q96AJ1
Cytogenetics:	16p13.3
MW:	48.2 kDa



[View online »](#)

Gene Summary:

The protein encoded by this gene contains a single coiled-coil region. Alternative splicing results in multiple transcript variants and protein isoforms. [provided by RefSeq, Jul 2012]