

## Product datasheet for **RC213513L4V**

### TSPAN17 (NM\_001006616) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	TSPAN17 (NM_001006616) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TSPAN17
Synonyms:	FBX23; FBXO23; TM4SF17
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001006616
ORF Size:	1109 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213513).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001006616.1</a>
RefSeq Size:	2573 bp
RefSeq ORF:	792 bp
Locus ID:	26262
UniProt ID:	<a href="#">Q96FV3</a>
Cytogenetics:	5q35.2
Protein Families:	Transmembrane
MW:	29.7 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a member of the transmembrane 4 superfamily. It is characterized by four tetraspanin transmembrane segments. The function of this gene has not yet been determined. [provided by RefSeq, Mar 2014]