

Product datasheet for **RC213520L2V**

KLHL15 (NM_030624) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	KLHL15 (NM_030624) Human Tagged ORF Clone Lentiviral Particle
Symbol:	KLHL15
Synonyms:	HEL-S-305
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_030624
ORF Size:	1813 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213520).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_030624.1
RefSeq Size:	6290 bp
RefSeq ORF:	1815 bp
Locus ID:	80311
UniProt ID:	Q96M94
Cytogenetics:	Xp22.11
MW:	69.8 kDa



[View online »](#)

Gene Summary:

This gene encodes a member of the kelch-like family of proteins that share a common domain structure consisting of an N-terminal broad-complex, tramtrack, bric-a-brac/poxvirus and zinc finger domain and C-terminal kelch repeat motifs. The encoded protein may be involved in protein ubiquitination and cytoskeletal organization. [provided by RefSeq, Apr 2009]