

Product datasheet for **RC213567L3V**

SLC14A2 (NM_007163) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | SLC14A2 (NM_007163) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | SLC14A2 |
| Synonyms: | hUT-A6; HUT2; UT-A2; UT2; UTA; UTR |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_007163 |
| ORF Size: | 2760 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC213567). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_007163.2 , NP_009094.2 |
| RefSeq Size: | 4166 bp |
| RefSeq ORF: | 2763 bp |
| Locus ID: | 8170 |
| UniProt ID: | Q15849 |
| Cytogenetics: | 18q12.3 |
| Domains: | UT |
| Protein Families: | Transmembrane |



[View online »](#)

MW: 101 kDa

Gene Summary: The protein encoded by this gene belongs to the urea transporter family. In mammalian cells, urea is the chief end product of nitrogen catabolism, and plays an important role in the urinary concentration mechanism. This protein is expressed in the inner medulla of the kidney, and mediates rapid transepithelial urea transport across the inner medullary collecting duct. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Sep 2011]