

Product datasheet for RC213569A1V

Human ASB13 (NM_024701) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human ASB13 (NM_024701) AAV Particle
Tag: Myc-DDK
Symbol: ASB13
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC213569 representing NM_024701
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGAGCCCCGGGCGGCGGACGGCTGCTTCCTGGGCGACGTGGGTTTCTGGGTGGAGCGGACCCCTGTGC
 ACGAGGCAGCCCAGCGGGGTGAGAGCCTGCAGCTGCAACAGCTGATCGAGAGCGGCGCCTGCGTGAACCA
 GGTACCGTGGACTCCATCACGCCCCTGCACGCAGCCAGTCTGCAGGGCCAGGCGCGGTGTGTGCAGCTG
 CTGCTGGCGGCTGGGGCCAGGTGGATGCTCGAACATCGACGGCAGCACCCCGCTCTGCGATGCCTGCG
 CCTCGGGCAGCATCGAGTGTGTGAAGCTCTTGCTGTCTACGGGGCCAAGGTCAACCCTCCCCTGTACAC
 AGCGTCCCCCTGCACGAGGCCCTGCATGAGCGGAGTTCGAATGTGTGAGGCTTCTATTGACGTGCGG
 GCCAATCTGGAAGCGCACGATTGCCATTTGGGACCCCTCTGCACGTTGCCTGTGCCCGGAGCATCTGG
 ACTGTGTCAAAGTGCTGCTCAATGCAGGGGCCAACGTGAATGCGGCAAAGCTTCATGAGACTGCCCTTCA
 CCACGCGGCCAAGGTCAAGATGTTGACCTCATCGAGATGCTTATCGAGTTTGGCGGCAACATCTACGCC
 CGGGACAACCGCGGAAGAAGCCGTCTGACTACACGTGGAGCAGCAGCGCTCCCGCCAAGTGCTTCGAGT
 ACTACGAAAAGACACCTCTGACTCTGTACAGCTCTGCAGGTGAACTTGAGGAAGGCCACTGGCGTCCG
 AGGGCTGGAGAAGATTGCCAAGTTAAACATCCCGCCCCGGCTCATTGATTACCTCTCTACAAC

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence:	>RC213569 representing NM_024701 Red=Cloning site Green=Tags(s) MEPRAADGCFLGDVGFWVERTPVHEAAQRGESLQLQQLIESGACVNQVTVDSTITPLHAASLQGQARCVQL LLAAGAQVDARNIDGSTPLCDACASGSIECVKLLL SYGAKVNPPLYTASPLHEACMSGSECVRLIDVG ANLEAHDCHFGTPLHVACAREHLDCKVLLNAGANVNAAKLHETALHHAAKVKNVDLIEMLIEFGGNIYA RDNRGKKPSDYTWSSAPAKCFEYYEKTPLTSLQCRVNLRKATGVRGLEKIAKLNIPRLIDYLSYN TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_024701
ORF Size:	834 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_024701.2
RefSeq Size:	2764 bp
RefSeq ORF:	837 bp
Locus ID:	79754
UniProt ID:	Q8W XK3
Cytogenetics:	10p15.1
MW:	29.8 kDa