

Product datasheet for **RC213610**

CHM (NM_000390) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CHM (NM_000390) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	CHM
Synonyms:	DXS540; GGTA; HSD-32; REP-1; TCD
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RC213610 representing NM_000390
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGCGGATACTCTCCCTTCGGAGTTTGATGTGATCGTAATAGGGACGGGTTTGCCTGAATCCATCATTG
CAGCTGCATGTTCAAGAAGTGGCCGGAGAGTTCTGCATGTTGATTCAAGAAGCTACTATGGAGGAACTG
GGCCAGTTTTAGCTTTTCAGGACTATTGCTGCTAAAGGAATACCAGGAAAACAGTGACATTGTAAGT
GACAGTCCAGTGTGGCAAGACCAGATCCTGAAAATGAAGAAGCCATTGCTCTTAGCAGGAAGGACAAAA
CTATTCACATGTGGAAGTATTTGTTATGCCAGTCAGGATTTGCATGAAGATGTCGAAGAAGCTGGTGC
ACTGCAGAAAAATCATGCTCTTGTGACATCTGCAAACTCCACAGAAGCTGCAGATTCTGCCTTCTGCCT
ACGGAGGATGAGTCATTAAGCACTATGAGCTGTGAAATGCTCACAGAACAACTCCAAGCAGCGATCCAG
AGAATGCGCTAGAAGTAAATGGTCTGAAGTGACAGGGGAAAAAGAAAACCATTGTGATGATAAACTTG
TGTGCCATCAACTTCAGCAGAAGACATGAGTGAATGTGCTATAGCAGAAGATACCACAGAGCAACCA
AAGAAAAACAGAATTACTTACTCACAATTTATTAAGAAGGCAGGAGATTAATATTGATTTAGTATCAA
AGCTGCTGATTCTCGAGGATTACTAATTGATCTTCTAATCAAATCTAATGTTAGTCGATATGCAGAGTT
TAAAAATATTACCAGGATTTCTGCATTTTCGAGAAGGACGAGTGGAAACAGGTTCCGTGTTCCAGAGCAGAT
GTCTTTAATAGCAACAACCTTACTATGGTAGAAAAGCGAATGCTAATGAAATTTCTTACATTTTGTATGG
AATATGAGAAATATCCTGATGAATATAAAGGATATGAAGAGATCACATTTTATGAATATTTAAAGACTCA
AAAATTAACCCCAACCTCCAATATATTGTCATGCATTCAATTGCAATGACATCAGAGACAGCCAGCAGC
ACCATAGATGGTCTCAAAGCTACCAAAAACCTTCTTCACTGTCTTGGGCGGTATGGCAACACTCCATTTT
TGTTTCCCTTATATGGCCAAGGAGAACCTCCCCAGTGTCTTGCAGGATGTGTGCTGTGTTTGGTGGAA
TTATTGCTTCGCCATTCAGTACAGTGCCTGTAGTGGACAAAGAATCCAGAAAATGTAAGCAATTATA
GATCAGTTTGGTCAGAGAATAATCTCTGAGCATTCTCCTGCGGAGGACAGTTACTTTCCTGAGAATGT
GCTCACGTGTGCAATACAGGCAGATCTCCAGGGCAGTGTGATTACAGATAGATCTGCCTAAAAACAGA
TTCAGATCAACAGATTTCCATTTTGCAGTGCAGCAGAGGAACCGAACTTTTGTGCTCGGGTCATT
GAGTTATGTTCTTCAACGATGACATGCATGAAAGGCACCTATTTGGTTCATTTGACTTGCACATCTTCTA
AAACAGCAAGAGAAGATTTAGAATCAGTTGTGCAGAAATGTTTGTTCATATACTGAAATGGAGATAGA
AAATGAACAAGTAGAAAAGCCAAGAATTCTGTGGGCTCTTACTTCAATATGAGAGATTCGTGAGACATC
AGCAGGAGCTGTTATAATGATTTACCATCAACGTTTATGTCTGCTCTGCCCAGATTGTGGTTTAGGAA
ATGATAATGCAGTCAAACAGGCTGAAACACTTTTCCAGGAAATCTGCCCAATGAAGATTTCTGTCCCC
TCCACCAATCCTGAAGACATTATCCTTGATGGAGACAGTTTACAGCCAGAGGCTTCAGAATCCAGTGCC
ATACCAGAGGCTAACTCGGAGACTTTCAAGGAAAGCACAAACCTTGAAACCTAGAGGAGTCTCTGAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC213610 representing NM_000390
Red=Cloning site Green=Tags(s)

MADTLPSEFDVIVIGTGLPESIIAAACSRSGRRVLHVDSRSYYGGNWSFSGLLSWLKEYQENS DIVS
 DSPVWQDQILENEEAIALSRKDKTIQHVEVFCYASQDLHEDVEEAGALQKNHALVTSANSTEADSAFLP
 TEDESLSTMSCEMLTEQTPSSDPENALEVNGAEVTGEKENHCDDKTCVPSTSAEDMSENVPIAEDTTEQP
 KKNRITYSQIIEGRRFNIDLVSKLLYSRGLLIDLLIKSNVSRYAEFKNI TRILAFREGRVEQVPCSRAD
 VFNSKQLTMVEKRMLMKFLTFCMEYEKYPDEYKGYEITFYEYLKTQKLPNLQYIVMHSI AMTSETASS
 TIDGLKATKNFLHCLGRYGNTPFLFPLYGQGELPQCFRCMCAVFGGIYCLRHSVQCLVVDKESRKCKAII
 DQFGQRIISEHFLVEDSYFPENMCSRVRQYRQISRAVLITDRSVLKTDSDDQIISILTVPAEEP GTF AVRVI
 ELCSSMTCKMGTYL VHLTCTSSKTAREDES VVQKLFVPYTEMIEENEQVEKPRILWALYFNMRDSSDI
 SRSCYNDLPSNVVYVCSGPD CGLGNDNAVKQAETL FQEICPNEDFCPPPPNPEDIILDGDSLQPEASSESA
 IPEANSETFKESTNLGNLEESSE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_000390

ORF Size: 1959 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_000390.4
RefSeq Size:	5442 bp
RefSeq ORF:	1962 bp
Locus ID:	1121
UniProt ID:	P24386
Cytogenetics:	Xq21.2
Domains:	GDI
Protein Families:	Druggable Genome
MW:	73.5 kDa
Gene Summary:	<p>This gene encodes component A of the RAB geranylgeranyl transferase holoenzyme. In the dimeric holoenzyme, this subunit binds unprenylated Rab GTPases and then presents them to the catalytic Rab GGTase subunit for the geranylgeranyl transfer reaction. Rab GTPases need to be geranylgeranylated on either one or two cysteine residues in their C-terminus to localize to the correct intracellular membrane. Mutations in this gene are a cause of choroideremia; also known as tapetochoroidal dystrophy (TCD). This X-linked disease is characterized by progressive dystrophy of the choroid, retinal pigment epithelium and retina. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Mar 2016]</p>

Product images:



Circular map for RC213610