

Product datasheet for **RC213633L3V**

TMEM98 (NM_001033504) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TMEM98 (NM_001033504) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TMEM98
Synonyms:	TADA1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001033504
ORF Size:	678 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213633).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001033504.1 , NP_001028676.1
RefSeq Size:	1732 bp
RefSeq ORF:	681 bp
Locus ID:	26022
UniProt ID:	Q9Y2Y6
Cytogenetics:	17q11.2
Protein Families:	Transmembrane
MW:	24.6 kDa



[View online »](#)

Gene Summary:

This gene encodes a transmembrane protein. A missense mutation in this gene result in Nanophthalmos 4 (NNO4). Alternatively spliced transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Sep 2014]