

## Product datasheet for **RC213683L1V**

### Cell adhesion molecule 4 (CADM4) (NM\_145296) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Cell adhesion molecule 4 (CADM4) (NM_145296) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Cell adhesion molecule 4
Synonyms:	IGSF4C; Necl-4; NECL4; synCAM4; TSLL2
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_145296
ORF Size:	1164 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213683).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_145296.1</a>
RefSeq Size:	2176 bp
RefSeq ORF:	1167 bp
Locus ID:	199731
UniProt ID:	<a href="#">Q8NFZ8</a>
Cytogenetics:	19q13.31
Protein Families:	Transmembrane
MW:	42.6 kDa



[View online »](#)

**Gene Summary:**

Involved in the cell-cell adhesion. Has calcium- and magnesium-independent cell-cell adhesion activity. May have tumor-suppressor activity.[UniProtKB/Swiss-Prot Function]