

Product datasheet for **RC213703L3V**

PCDH1 (NM_032420) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PCDH1 (NM_032420) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PCDH1
Synonyms:	PC42; PCDH42
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_032420
ORF Size:	3711 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213703).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_032420.2
RefSeq Size:	4779 bp
RefSeq ORF:	3714 bp
Locus ID:	5097
UniProt ID:	Q08174
Cytogenetics:	5q31.3
Domains:	CA
Protein Families:	Transcription Factors, Transmembrane



[View online »](#)

MW: 130.1 kDa

Gene Summary: This gene belongs to the protocadherin subfamily within the cadherin superfamily. The encoded protein is a membrane protein found at cell-cell boundaries. It is involved in neural cell adhesion, suggesting a possible role in neuronal development. The protein includes an extracellular region, containing 7 cadherin-like domains, a transmembrane region and a C-terminal cytoplasmic region. Cells expressing the protein showed cell aggregation activity. Alternative splicing occurs in this gene. [provided by RefSeq, Jul 2008]